

Shexie

Overview

This guide provides instruction on configuring both Medical-Objects Capricorn and Shexie software to import electronic results.

If you require any assistance with the installation or configuration, or if you have any questions, please contact the Medical-Objects Helpdesk on (07) 5456 6000.

If you have not installed the Medical-Objects Download client, please see the guide for [Installing Capricorn on Windows](#).

You can also contact the Shexie helpdesk on 1300 SHEXIE (1300 743 943)

- 1 Overview
- 2 Configuring Capricorn Software
 - 2.1 Accessing the Capricorn Configuration Window
 - 2.2 Capricorn Configuration Settings
- 3 Shexie Configuration

Configuring Capricorn Software

Accessing the Capricorn Configuration Window

Note:

In some cases, the Capricorn may be configured as a service. See [this guide](#) for instructions on how to launch the Capricorn if it is set up as a service.

1. To launch the HL7 Tray application, Go to **Start Menu, All Programs, Medical-Objects** and select the **Medical-Objects Capricorn**. Alternatively, You can go to **C:\MO\Capricorn** and run the **Capricorn.exe**.
2. The following icon will appear in the tray. Double click on it to bring up the log window.




3. When the log window is displayed, Click on the **Configuration** icon highlighted below, or select **Utility** and select **Configuration**.



Capricorn Configuration Settings

1. Click on the **Integration Options**.



2. Set the **ACK Management** path to the local Import folder on the server or the network path on a workstation, that points to the server import folder.
For example:

Acknowledgement paths in Shexie

```
(Capricorn installed on the Server) --> C:\Shexie\pathology\MoAcks
Or
(Capricorn installed on the Workstation) --
> \\server\Shexie\pathology\MoAcks
```

Note:


If the folder does not exist, you must create it.

Integration Settings

Set Tick Box to **"Evaluate and Strip"**

Set Incoming Modifier to **"Shexie Modifier"**

3. Select **Server Parameters** from the left panel.



4. Under **Message Output Directory (main)** set the results import path.
5. The path can be different depending on where Capricorn is installed.

For example:

Download paths in Shexie

```
(Capricorn installed on the Server) --> C:\Shexie\pathology\  
Or  
(Capricorn installed on the Workstation) --  
> \\server\Shexie\pathology\
```


6. The configuration is complete. Click **Apply** and **OK**.

Shexie Configuration


Note:

If performing the steps below doesn't show any of your reports and they stay remaining in the folder - Call Shexie (1300 743 943) as they will need to further configure Shexie to Import.


1. To check Pathology results, go to the **Housekeeping** Tab and then select **Pathology Results**.




2. This will then display the Pathology Results you have received.



3. To check if Specialist or Allied health letters have arrived select **E-Letters**



4. This will then display the E-Letters you have received.



5. You can also Manually force import results, which can be done from the Windows taskbar.

