

# Relay Client - Eclipse

- [System Requirements](#)
- [How to setup a basic send queue](#)
- [Management and Monitoring tasks for Eclipse](#)
- [Saving Application Acks for your PMS or RIS](#)
- [Extra Scripts](#)
- [Interfacing with the Sharp Scanner](#)
- [Installing Eclipse on Windows](#)

## Overview

Eclipse is Medical-Objects' integration product.

Based on a modular design it allows for manipulation of health messages so that many actions can be taken on a message so that it will be the destination format.

It is useful for integration to other clinical tools such as **Practice Management Systems** and Clinical Systems such as **Cerner** and **IBA**.

Eclipse allows you to do a number of things, including but not limited to:

- Convert a HL7 message to PIT
- Convert a PIT message to HL7
- Convert a HL7 message to Text
- Email the output of the message

Eclipse includes approximately 100 modules that you can use to modify the message.