

# How to access Capricorn Configuration

## Overview


This guide provides basic instruction on how to find and launch the Medical-Objects Download Client (Capricorn) and access the configuration panel.

- 1 [Overview](#)
- 2 [Launching the Capricorn From Services](#)
- 3 [Launching the Capricorn, and the Configuration Window](#)


## Launching the Capricorn From Services

If the Medical Objects Download Client (Capricorn) icon is not in the tray, It may be running as a service.


1. Go to your Start menu, and type **Run**, and Select the **Run** application.
2. In the box that appears, Type **Services.msc** and click **OK**. This will open the Services window.



3. Look for the **CAPRICORN** service on the list. If it is not on this list, Capricorn is not installed as a service. In this case, Skip to the next section.



4. Right click on the **CAPRICORN** service, and select the **Stop** option.



## Launching the Capricorn, and the Configuration Window

1. To launch the HL7 Tray application, Go to **Start Menu, All Programs, Medical-Objects** and select the **Medical-Objects Capricorn**. Alternatively, You can go to **C:\MO\Capricorn** and run the **Capricorn.exe**.
2. The following icon will appear in the tray. Double click on it to bring up the log window.



3. When the log window is displayed, Click on the **Configuration** icon highlighted below, or select **Utility** and select **Configuration**.

