Direct CONTROL

Overview

This article is to assist with the setup of Medical Objects and Direct Control

Requirements

- ①1. Direct CONTROL version 22.2.7.1 or higher is required. This can be checked in the top right corner of Direct CONTROL
 - You are in Direct CONTROL. Version: 22.2.7.1
 - 2. Capricorn 3.2.13.138 or higher is required.

Configuring Capricorn Software

Accessing the Capricorn Configuration Window

In some cases, the Capricorn may be configured as a service. See this guide for instructions on how to launch the Capricorn if it is set up as a service.

- To launch the HL7 Tray application, Go to Start Menu, All Programs, Medical-Objects and select the Medical-Objects Capricorn. Alternatively, You can go to C:\MO\Capricorn and run the Capricorn.exe.
- 2. The following icon will appear in the tray. Double click on it to bring up the log window.



3. When the log window is displayed, Click on the **Configuration** icon highlighted below, or select **Utility** and select **Configuration**.



Capricorn Configuration Settings

 Click on Ack Management. Set the Folder for Acks created by this practices PMS (Receiver) to the HL7 Out path in the walkthrough here: HL7Out 1 Overview 2 Requirements 3 Configuring Capricorn Software 3.1 Accessing the Capricorn Configuration Window 3.2 Capricorn Configuration Settings 4 Direct CONTROL Import Configuration 4.1 HL7 In 4.2 HL7 Out 5 Inbox



2. Click on Integration Options.

3.1 Check that the Incoming Modifier is set to Direct CONTROL



3. Select Server Parameters from the left panel.



The paths above are set as follows:



If the folder does not exist, You must create it.

4. The Medical-Objects configuration is complete. Click Apply and OK.

Direct CONTROL Import Configuration

In Direct CONTROL navigate to Utilities Configure

HL7 In



If paths are not populated click "Load Defaults" then "Install/Upgrade"

Change settings as required.

Set the download path in HL7 In Location as the first setting in HL7 Basic Setup.

Select "Start" button

HL7 Out

This is currently only for return acknowledgements sending messages directly from Direct CONTROL has not been finalised.

Configuration										_
Advanced Digital Services File Loc	ations General HL7	Installation Integrations	International	Medicare	Menus	Portal	Printing/Scanning	SMS	٠	ি
HL7 In HL7 Out										Suppor
	General Configuration									
Save HL7 files to:										
	MSH Configuration									
Receiving Application										
Sending Facility:										
	Acknowledgement Mess	lages								
Save Acknowledgement messages to:	E:\DCDEMO4\DirectContr	olServerFiles\HL7\HL7Out\								

Set the "Save Acknowledgement messages to:" path.

Inbox

Regular Tasks --> Communications/Results

If Direct CONTROL was not able to find a match for an existing Client, the documents will appear in the Unabcated tab.

Communication Management																
	Entity:	(AI)			~	Status:	(AI)			lode 🖌	1					
Delivery	Method:	(AI)				¥										
Fro	om Date:	Tuesday, 6 September 2022				•	To Date:	: 🔽 Thursday, 6 October 2022 🔹					Load Data			
All) Administration Admissions Calls Care Plans Clinical Custom Emails Messages Referrals Results Unaflocated SMS Transactions [Delet																
ClientName Communi			ation Name					DateEntered		DateOfCommun	ication	Comments	ci 🔺			
Co	ntains:		7	Contains:					Y	Equals:	7	Equals:	7	Contains: 🌱		
+ Mo	use, Mick	ey		LETTER -	DISCHARGE	SUMMAR	RY-pms3me	dd.PDF		31/08/2022 8:48	MA 00:	26/08/2022 4:41	:00 PM		-	
Mo	ouse, Mick	ey		LETTER -	Corrected DI	SCHARGE	E SUMMAR'	Y-pms3me	dd.PDF	31/08/2022 8:48	MA 00:	26/08/2022 4:41	:00 PM			