

Installing the Referral Client

Overview

This guide provides an overview of the installation process for the medical-objects Referral Client (or Trinity). The software can also be configured with your PKI token. Please see [this guide](#) after installation is completed.


- 1 [Overview](#)
- 2 [Installation and configuration of Medical Objects Trinity](#)
- 3 [Large deployments](#)

Installation and configuration of Medical Objects Trinity


- 1. Download the Referral Client Setup from the following link.

<https://download.medical-objects.com.au/>

- 2. The link on the download page looks like the following.



- 3. The following window will appear.



- 4. Click on **Save**.

- 5. The **Save As** window will appear. Save the **referral-setup.exe** in the directory **C:\MO**.

Note:


 If the directory does not exist, Create it.

- 6. Once you have clicked **Save**, the software will download.




- 7. After it has been downloaded click **Run**. If a security warning comes up, Click **Run** again.


- 8. Once you have clicked **Run**, the installation process will begin.




9. Click **Next**.



10. Click **Next**.



11. Click **Install**.



12. Click **Finish**.

The setup is now complete. The application should automatically launch in a few seconds.

Note:



For Local networks, the following firewall rules will need to be put in place.

UDP Port 21734 OUT on the Referral Client (Trinity) PC

UDP Port 21734 IN on the Download Client (Capricorn) PC

Large deployments

Note the switches `/SILENT` or `/VERYSILENT`, in conjunction with the installer, will suppress the installer dialogues partially or completely.