

# Upgrading the Medical-Objects Download Client - Mac

## Overview


This guide provides an overview on updating the Medical-Objects Download Client (Capricorn) For Mac based operating systems.

If you have any questions or require any further assistance, Please contact the Medical-Objects Helpdesk on (07) 5456 6000.

- 1 [Overview](#)
- 2 [Checking your version](#)
- 3 [Upgrading the software](#)

## Checking your version

The Medical-Objects Download Client can be upgraded by installing the newer application of the top of the existing installation. You can check your current version by checking the Activity report, as explained in this guide [here](#). The Capricorn version number will be listed as in the screenshot below.



## Upgrading the software

Note:




Make sure you have a backup of your software before you continue. If you are not sure how to make a backup, please see [this guide](#).

Before you upgrade the application, You will first need to make sure that it is no longer running.


1. Open the Medical-Objects Download Client by clicking the globe icon shown below located in your tray.




2. This will bring up the below window.



3. Click on the **Exit** button to shut down the Medical-Objects Download client.
4. You will then need to download the latest version from our website. Open the Web browser and navigate to [download.medical-objects.com.au](http://download.medical-objects.com.au)
5. Under Latest Installers, Find **Download Client (Capricorn)** and Click on '**Download (Mac)**'
6. From the drop down select '**OS X 10.4 and Higher**' to start the download of '**MO-Capricorn-Setup.dmg**'



7. Click the '**Downloads**' icon in safari (or navigate to the downloads folder)



8. Double click the '**MO-Capricorn-Setup.dmg**' in the downloads list (or from the downloads folder).
9. The Mac will verify the integrity of the DMG and then open a finder window displaying the **MO-Capricorn-Setup.app** file.
10. Double click on the '**MO-Capricorn-Setup.app**' to begin the installation.

Important:



After running the **MO-Capricorn-Setup.app** you may be confronted with a popup informing you that the **MO-Capricorn-Setup.app** has not been downloaded from the App store and only provide you with a '**Cancel**' button. if you get the popup, follow steps 10a-10c of this guide. If you don't get a popup and the Capricorns Token windows appears, Go to Step 11.


10a. If you get the pop-up saying 'The Capricorn has not been downloaded from the App Store', close the pop up and **Right-Click** on the **MO-Capricorn-Setup.app**.

10b. From the menu, Select **Open**.

10c. The popup box will appear again, but this time you can select the **Open** option.

Info:

The Pop-up described above is part of Apples 'Gatekeeper' security.



11. The Application should Automatically upgrade. When it is complete, the Download client will reopen.  
If you are unsure if the upgrade was successful, please contact the Medical-Objects helpdesk before you make a backup of the newly upgraded application.