

Overview

1 Overview

1. HPI - **History of Present Illness**

HPI - I

HPI - O

Pathology Requests	Radiology Requests	Other Requests	Referral/Documents Sent	Radiology Reporting
<input type="checkbox"/> Allow Pathology/Requests Request Form <input type="text"/> <input type="checkbox"/> Don't print requesting doctor details				
List order <input type="text"/>				
Test List <input type="text"/>				

[illegible]

1. Create the contact as normal, but leave the **Request form (Pathreqo)** and **list order** fields blank for now. This should let you create the base of the contact without the error message.
2. **Save, Exit** and **re-enter** the contact. You should now be able to setup the '**messages in**' part, HOWEVER you will find that when you are trying to select **ORU^RO1** as the message type, that it doesn't exist. If this is the case you should choose **ADT^ZHI Healthconnect Discharge Summary** as a temporary measure (see Below) Now save the messages in part and exit the contact completely

Receive Messages From

Close

Receive From: (Select)

Description:

Inactive: ☐

Sending Application (or PIT Source):

Sending Facility:

Transport Type: (None)

TCP Port Number:

Protocol:

Folder Name:

Filename Mask:

Acknowledgement Type:

Ack Folder Name:

Ack Filename Scheme: (None)

Message Type

Select the Message Type

(Select)

(Select)

XQ-3241 - HealthConnect Discharge Summary

REF2743 - HealthConnect Discharge Summary

REF4112 - Patient Referral/Genetic Report

Message Details

Message Type	Message Type Description	Check Run No.

Add

Remove

3. Go Back into the contact and you should now be able to input the **Request form (Pathreqo)** and **list order** fields without an error. **Save** and exit the contact.
4. Re-enter the contact and go back into the '**messages in**' section. You should now be able to select **ORU^RO1** as the message type. Delete the **ADT^ZHI** part that you chose before.
5. save and exit. The contact creation is now complete.

