

# Pathology Printing Setup

Templates for Pathology Printing can be found in a standard Explorer install at C:\MO\Clients\CSS\_TEMPLATES

1 [Fields for Reference](#)

## Fields for Reference

DOB

PatientName

RequestedBy

Requested

TestDate

Reported

Printed

PatientHomePhone

PatientHomeMobile

PatientHomeFax

PatientAddressFull

PatientAddress

PatientAddressWithHTMLBreak

PatientAge

PatientSex

YourRef

LabNo

CopiesToHorizontal

CopiesToHorizontalWithoutCurrentXCN

SelectedXCNListWithoutCurrentXCN

CopiesToVerticle

StyleSheetMOPrint

StyleSheetChangelayout

TargetDoctorFullAddress

TargetDoctorFullAddressCapCase

TargetDoctorFullAddressWithProviderID

TargetDoctorsCCList

PageNumber

PageCount

CopyFor

LABDISPLAY

HIGHLIGHTCOLOUR

REPORTEDBY

CurrentXCN

URNUMBER

OriginalRecipient