

Marking a Requested Chart as Reviewed

Overview

- 1 Overview
- 2 Marking a Requested Chart as Reviewed

This guide provides a brief overview on marking a requested charged as reviewed within the Medical-Objects Explorer Client.

If you require any assistance with the installation or configuration, or if you have any questions, please contact the Medical-Objects Helpdesk on (07) 5456 6000.

Marking a Requested Chart as Reviewed

To mark a Requested Chart as Reviewed, click on the Requested Charts Folder to display the list of Patient's Charts that have been requested by Doctors within the Practice.

Patients

Query

Unreviewed Results

Further Review

Requested Charts

SMITH, JOHN

Decision Support

Reports created

0

1

1

1

0

0

Medical Objects - Logged into Test Message application as THOMAS WEBER (THW45500162)

Patients

Available Results

GP Referral

John SMITH

22 Imagination Lane 8015BAYVIEW QLD 4001

GP Referral

Medical-Objects Test Surgery

Suite 10 / 180 Sussex Road

Norrockyshore, QLD, 4358

Tuesday, 12 May 2014

Dr Jane DOE (J4545888202)

Dr Jane DOE (J4545888202)

Age 35years Medical Clinic

123 Test Street

Brinsford, QLD, 4000

AL: Mr John Smith

DOB: 15/11/1981

22 Imagination Lane

Bayview, QLD 4001

Dear Dr Doe,

Thank you for referring Mr John Smith for our services. We have found that Mr Smith

has no serious illness, however he could benefit greatly from eating more healthy

foods.

Mr Patient was concerned as to what foods to eat, so I was able to confirm for him

the following foods:

1 Vegetables

2 Fruits

Cancel Chart

Once a Patient's Chart has been flagged as a Requested Chart, the Cancel Chart Request option appears on the Patient's Results screen. Use this option to cancel the Chart request, or acknowledge the Chart has been provided.