

Viewing Results Marked for Further Review

Overview

- 1 [Overview](#)
- 2 [Viewing Results marked for Further Review](#)

This guide provides a brief overview on Viewing Results Marked for Further Review using the Medical-Objects Explorer (Meridian).

If you require any assistance with the installation or configuration, or if you have any questions, please contact the Medical-Objects Helpdesk on (07) 5456 6000.

Viewing Results marked for Further Review

To view Results marked for Further Review, Click the Further Review Folder in the Patients Box on the left hand side of the screen.

Patients	
Query	0
Unreviewed Results	1
Further Review	1
SMITH, JOHN	
Requested Charts	0
Decision Support	0
Reports created	0

The folder displayed shows the number of Patients that have Results marked for further Review. To Review a patients Results, Click on the Patients Name and the number of Results that have been Parked for Further review will appear in the Available Results box on the left hand side of the screen. The Result that has been parked will appear on the right hand side of the screen. For more information on Parked Results, See [This guide](#).

Available Results	
All	Disease
Groups	Graphs
21.05.2014 : GP Referral	