# **Provider Mail Merge**

### Overview

The Provider Mail Merge functionality in Explorer allows users to perform printing tasks specific to providers. It can be accessed from Provider List For Mail Merge under the Providers menu.

File Patients Documents Create Document	Providers View Navigate Tools Reports Window Help
🕼 🔕 📄 💿 🧇 🔊	P Lookup Provider Provider Map
	Edit my Profile
	Providers online Advanced Provider List Provider List Provider List For Mail Merge

1 Overview 2 Selecting Providers 3 Output Types 3.1 Print List 3.2 Print Labels 3.3 Merge Print

Here you are able to search for and select a list of desired providers and then print that list in three different ways, including as a list, as labels and as mail merged letters.

## **Selecting Providers**

First, search for providers by entering either one postcode or multiple postcodes separated with commas.

You may also drop down the box to select previous postcodes that have been entered.

After entering one or more postcodes, click the Lightning Bolt button to begin the search.

A list of providers within the given postcode(s) will be returned, detailing their name, practice name, address and provider ID.

Use the check boxes next to the providers' names to select/deselect the providers you want to continue with.

Enter postcode(s	) here Beg	Search Provider Mail Merge		×
Postcode(s) ( either 4000 or	4001,4002,4003)	Filter by Job Class		
4558	<ul> <li>Within radius (km)</li> <li>30 C</li> </ul>			
Nare	Practice Name 3ob	lass Addre i	Provider ID	^
	Extend or minimise	Select between GP	-	
	search from postcode(s)	Specialist or Staff		
REPORTED THEY	THE REACTION	TOT IN ACTIVE MAN COLUMN DE 4559		
Select the Provid	ters you wish to send to 1		_	
		The latter was an end of the second sec		
¥		No. of a local distance of the local distanc		
2				
		The second se		
12 Desiders Frank II Column				× 1
Literiondes round (3 seecles				
outpu	t types	Prevenu		
Merce From Texniste	Phil Outpu			
Specialist Letter	<ul> <li>Setup</li> </ul>	There are a second seco		
Print to		Proview		
	<ul> <li>Setup</li> </ul>			
	Print			
			Close	

## **Output Types**

There are three formats you can print the provider details in and they are as follows:

#### Print List

Print list prints a list of the selected providers - including their name, practice name, address and provider ID - in a table on a landscape-oriented page.

Use the drop down list to select the printer you want to print from and click the Setup button if you want to configure your printers.

Click the Print button to print the list.



#### **Print Labels**

Print labels allows you to print a sheet of labels that contain the providers' details (one label for each selected provider).

Select the label you want to use from the drop down box and click the Print button to print the labels.

Output Types	Preview
Print Lat Print Labels Marge Print	Copies: 1 🖉 Start at label: 1 💇 🔇 🔕 🔕 🔇 🕲 Page 1 of 1
kteri 3 • Setup	Open Label Configuration
Print Label	Dr. Medical Cantre
	BRISBANE OLD 4001

#### Note

Note: Only labels that are set up with ProviderAddress will show in the list. To setup an existing label or create a new label, click the Edit Label Config button.

<b>6</b>	Label Configuration	×
Selford Layout Hame Lobb 1 Prese (W. 283, H1002) Test Devention De	Single Label Preview	ostoa
Test Print		OK Cancel

#### Merge Print

Merge print allows you to print letters, based on a template, for each of the selected providers with the providers' details merging with the template's merge fields.

Select an appropriate template from the Merge From Template drop down list and then select the printer you want to use from the Print to list. Click the Print button to print the letters.

You may also open the template setup and printer setup windows by clicking the appropriate Setup buttons.

Note: Note: Provider name, practice, address and provider ID will only appear if the selected template has those specific merge fields in them.

