

# Reviewing Results

## Overview

- 1 [Overview](#)
- 2 [Reviewing results](#)


This guide provides an overview on reviewing Results within the Medical-Objects Explorer (Meridian).

If you require any assistance with the installation or configuration, or if you have any questions, please contact the Medical-Objects Helpdesk on (07) 5456 6000.

## Reviewing results

There are several ways to review results using Explorer.

1. To Review Results, go to the **Patients** box at the top left hand side of the screen.



Query	0
Unreviewed Results	1
Further Review	0
Requested Charts	0
Decision Support	0
Reports created	0

Folder Meanings - If you click on a folder once it will expand or collapse.

**Query** - Shows the number of queries regarding Patients that have been performed, and when the Query Folder is expanded it will list the Results of each **query**.


**Unreviewed Results** - Shows the number of Results that have been sent to you that are still not marked as **reviewed**.

**Further Review** - Displays the number of Results that have been "parked" awaiting further **review**. For more information on this see [this link](#).

**Requested Charts** - Shows the number of Patient's Charts that have been requested by any Doctor in the **Practice**.

**Decision Support** - Information that could help you make Clinical Decisions eg- keep information on a specific drug from MIMS on hand.


2. New Results are automatically downloaded by Explorer. However, Results can be retrieved to ensure that there are no outstanding Results .They can be retrieved in one of two ways:
  - a. By selecting **Unseen Documents** from the **Tools** Menu
  - or
  - b. By simply clicking on the **Binoculars** button in the top left of the Toolbar.



- When there are New Results, the number of Patients with new Results will be displayed next to the **Unreviewed Results** folder in the Patients Window. To retrieve the Results, Click on the **Unreviewed Results** folder (once only) to show the List of Patients:

Query	0
Unreviewed Results	1
SMITH, JOHN	
Further Review	0
Requested Charts	0
Decision Support	0
Reports created	0

- Click on the individual Patient's Name to display the List of Results to be reviewed in the **Available Results** box.



- Click on an item in the list to display the Result. The lightning sign will appear next to the result to indicate that the Result is unreviewed. On the right hand side of the screen the result will be displayed. To easily scroll through the results, hit the space bar.

Tab Meanings.

