BreastScreen NSW Acknowledgements

Overview

This article is for customers of BreastScreen NSW and Best Practice. The focus of this article is configuring practice software result import acknowledgements to BreastScreen NSW so they know Best Practice has imported the messages.

Configuration

This is no longer needed if you have 1.11.1.934 and above of BP and have set up importing per the new setup instructions here.

Ack Path

Note:

In some cases, the Capricorn may be configured as a service. See this guide for instructions on how to launch the Capricorn if it is set up as a service.

- To launch the HL7 Tray application, Go to Start Menu, All Programs, Medical-Objects and select the Medical-Objects Capricorn. Alternatively, You can go to C:\MO\Capricorn and run the Capricorn.exe.
- 2. The following icon will appear in the tray. Double click on it to bring up the log window.



3. When the log window is displayed, Click on the **Configuration** icon highlighted below, or select **Utility** and select **Configuration**.



Capricorn Configuration Settings

1 Overview 2 Configuration 2.1 Ack Path 2.2 Capricorn Configuration Settings

1. Click on Ack Management.



2. Note the path set in ACK Management.

Note: If the folder does not exist, you must create it. Normally the default is C:\MO\Results\BP\Acks.

Make sure the results folder has **Modify, Read** and **Write** permissions, **AND** is shared on the Network.

3. Open the Best Practice software.



- 4. Click Setup from the toolbar shown above.
- 5. Click on Configuration from the drop down menu.



- 6. The **Configuration** window should be displayed.
- 7. Select Results Import from the left panel.
- 8. Select the Add button under Acknowledgements shown below.

, A	Configura	tion	
General General Results import Database	Configura Import investigation reports on this com Ont allow unallocated reports into holo Add header to incoming letters Report file search paths: IVPC2:VResults/XBPUIN Keep HL7 and PIT files after import Denta file selencing:	tion outer xx	Add Remove Check paths
Lists Reminders Care plans Care plans Email	Report the extension: CP CQP GR GR GR HI7 HMP HMP Momber of months to display in invastigatio Number of months to display in invastigatio Number of months to display in invastigatio Mumber of months to display in invostigatio Mumber of months to display in invastigatio Mumber of months to display in invastigatio	Accovedgements: Facility ACK p ACK p ACK p ist rup list: (G) reports list	Add Edit Edit Remove
	v		Save Cancel

- 9. The Send Acknowledgements window will appear.
- 10. Add BreastScreen NSW as the Facility name (include the space).

Send acknowledgeni	EIIIS	~
Facility name:	BREASTSCREEN NSW	
Acknowledgement path:	C:\MO\RESULTS\BP\ACKS	
	Same Car	aal

11. Click the file browse button to select the folder that was configured in Medical-Objects Capricorn in step 2. If you know what a UNC path is; then it is advised you use a UNC path versi on if you know it for the folder noted in step 2.

eg: \\PC3\Results\BP\Acks

Browse for Folder 🛛 💷	×
A Steverk	~
⊳ IN PC1	
⊿ [▲ PC2	
a 🎍 Results	
ACKS	
b Je Users	
⊳ 1 PC3	
▷ [♣ PC4	~
OK Cancel	

- 12. Once you have selected the path click **OK**.
- 13. Exit the configuration confirming that you have made changes. Best Practice should now return the acknowledgements that BreastScreen NSW require.