

# Capricorn Notification Settings

## Overview

- 1 [Overview](#)
- 2 [Configuration](#)

This guide contains the instructions to Configure the Medical-Objects Capricorn to show notifications when a new document is received.

If you have any questions or require any assistance, Feel free to contact the Medical-Objects Helpdesk on (07) 5456 6000.

## Configuration

Note:




In some cases, the Capricorn may be configured as a service. See [this guide](#) for instructions on how to launch the Capricorn if it is set up as a service.

1. To launch the HL7 Tray application, Go to **Start Menu, All Programs, Medical-Objects** and select the **Medical-Objects Capricorn**. Alternatively, You can go to **C:\MO\Capricorn** and run the **Capricorn.exe**.
2. The following icon will appear in the tray. Double click on it to bring up the log window.



3. In the Capricorn log window, tick the **Notify** Checkbox.



4. When a new document is received, a popup will display for several seconds as shown below.

