

About the Referral Client

Overview

The Medical-Objects Referral Client (Trinity) is a lightweight software program that can capture General Practitioner referrals from your **Practice Management Software** and send it to other connected Medical-Objects users.

The Referral Client is designed on standard network setups. While it can usually connect when using VPN's, the Medical-Objects Helpdesk does not provide support on non-standard network setups.


Diagram

Trinity broadcasts on UDP port 21734, if the broadcast is unobstructed the Capricorn will hear and respond with a UDP Response containing 'I'm Here' Details like IP and port number. Then communication commences via TCP.

Trouble Shooting

Set Firewall Application Rule for **Trinity.exe (Allow Outgoing all)**

Set Firewall Application Rule (**Allow Incoming all**)



- 1 [Overview](#)
- 2 [Diagram](#)
- 2.1 [Trouble Shooting](#)