


NASH

How to import via Capricorn:

- 1 [How to import via Capricorn:](#)
- 2 [How to import via the Referral Client \(Trinity\)](#)


1. Right click the Capricorn system tray and select **Configuration** from the context Menu.



Note:

If the Capricorn is not running in the system tray go to **Start** and select **Run**. the run box type in **services.msc**. Locate the service named '**Capricorn**' and stop the service. You can start the Capricorn as a tray application by going to **Start, All Programs, Medical-Objects, Capricorn**.

2. Select **PKI Configuration** and then hit the **Edit Keyring** button.




3. Click on the **Import Certificate File** button


If password issues persist, please contact Medicare.

How to import via the Referral Client (Trinity)


1. Right click the Trinity icon in the system tray.




2. Click on **Load PKI Location Certificate into Capricorn** from the **Admin** menu.
3. Click on the browse button for the signing certificate.



4. Browse to the **NASH Certificate CD** and import the **fac_sign** file.



5. Enter in the **PIC Passphrase** (password) for the NASH certificate and then hit **OK**.



The image shows a Windows-style dialog box titled "Load HESA certificates". It has a standard title bar with a minimize button, a maximize button, and a close button (red X). The dialog contains a section titled "Load PKI Certificates:" with three input fields. The first field is labeled "Encryption Certificate File" and is empty. The second field is labeled "Signing Certificate File" and contains the text "C:\Users\Wason\Desktop\fac_sign.p12". The third field is labeled "HIC Supplied Password" and contains a series of asterisks "*****" with a cursor at the end. This field is highlighted with a red rectangular box. At the bottom of the dialog, there are two buttons: "OK" and "Cancel". The "OK" button is also highlighted with a red rectangular box.

Load HESA certificates

Load PKI Certificates:

Encryption Certificate File

Signing Certificate File

C:\Users\Wason\Desktop\fac_sign.p12

HIC Supplied Password

OK Cancel