


# NASH

## How to import via Capricorn:


1. Right click the Capricorn system tray and select **Configuration** from the context Menu.



Note:

If the Capricorn is not running in the system tray go to **Start** and select **Run**. the run box type in **services.msc**. Locate the service named '**Capricorn**' and stop the service. You can start the Capricorn as a tray application by going to **Start, All Programs, Medical-Objects, Capricorn**.

2. Select **PKI Configuration** and then hit the **Edit Keyring** button.



3. Click on the **Import Certificate File** button


- 1 [How to import via Capricorn:](#)
- 2 [How to import via the Referral Client \(Trinity\)](#)




If password issues persist, please contact Medicare.

## How to import via the Referral Client (Trinity)


1. Right click the Trinity icon in the system tray.



2. Click on **Load PKI Location Certificate into Capricorn** from the **Admin** menu.
3. Click on the browse button for the signing certificate.



4. Browse to the **NASH Certificate CD** and import the **fac\_sign** file.



5. Enter in the **PIC Passphrase** (password) for the NASH certificate and then hit **OK**.

Load HESA certificates

Load PKI Certificates

Encryption Certificate File

Signing Certificate File

C:\Users\Wason\Desktop\fac\_sign.p12

HIC Supplied Password

OK Cancel